

# CAPOVARIO

**Our most versatile head protection helmet**

**Hernik**



Removable inner padding



Adjustable fit



Washable



## **Comfortable & secure**

We rely on the latest **technology** to make our products safe.

✓ Durable plastic structure for **maximum protection**



## **Trust in experience**

**20 years of experience** developing and manufacturing the best products.

✓ **Innovative mechanism** for head size adjustment



## **Made in Germany**

All our products are developed and **manufactured in Germany.**

✓ Certified according to **ISO 9001**





## Our CAPOVARIO

The CAPOVARIO is the most flexible of our head protection helmets. The anatomically shaped base structure made of solid plastic allows for head circumference adjustment (up to 4 cm).

The standard-integrated forehead protection shields the face and nose. If a flatter forehead protection is desired, the "LGT" version can be chosen as an alternative.


The helmet can be configured with various options tailored to it.

### Customize the helmet

The human head is as unique as a fingerprint. Therefore, great emphasis was placed on maximum flexibility during the development. The following can be adjusted with just a few simple steps:

- 1 Head Circumference
- 2 Carrying Straps
- 3 Crown Padding

Size	Head circumference	Weight
0	44 - 48 cm	270 g
1	48 - 51 cm	280 g
2	52 - 56 cm	290 g
3	56 - 60 cm	305 g
4	58 - 62 cm	320 g



*We are very proud of the many positive customer ratings. On our website you can find out, why Hernik regularly achieves top ratings on comparison portals.*

**Markus Jacksch,**  
Managing Director




## Custom accessories

To meet a wide range of requirements, innovative accessories were also considered in the development of the CAPOVARIO.

For example, we offer **ear protection (1)**, a **soft or hard closed top (2)**, or an **impact-resistant visor (3)**. Additionally, the helmet can be equipped with a modifiable **chin protector (4)** or a **mouth protection (5)**.

Additional accessories can be found in the corresponding order form.



## Material and colors

The CAPOVARIO uses a high-quality microfiber cover that is breathable, washable (40°C), and disinfectable, significantly improving hygiene.

The removable and washable inner padding also allows for comfortable air circulation thanks to its quilted design.

The helmet is currently available in 6 attractive colors.



We are certain that we can only develop the best head protection helmets together with our customers. That is why we are very happy to receive your suggestions for improvement and will incorporate them into our development.

Please send us a message:  
[office@hernik.de](mailto:office@hernik.de)





## Further products from HERNIK



### The right helmet for every need

Since every human head is unique, HERNIK offers various models to meet a wide range of needs.



### Child seats and seating shell concepts

Children with special needs require custom-designed restraint systems for safe transport in cars. We offer suitable car seats for children from toddler age up to adulthood.

We also provide custom seating shells, such as the aluminum shells ANATOMIC UNICO and ANATOMIC FLEX.

# Hernik

**Discover all our products  
on [hernik.de](http://hernik.de)**

**HERNIK GmbH**

Gewerbering 8,  
83543 Rott am Inn

+49 (08039) 9093930  
[office@hernik.de](mailto:office@hernik.de)  
[www.hernik.de](http://www.hernik.de)

