

# HERO

XXL

**Hernik**

**The rear-facing seat for big children**



~ 2-9 years

61 - 125 cm

up to 36 kg



## **Comfortable & secure**

We rely on the latest **technology** to make our products safe.

- ✓ **Best ratings** on comparison portals
- ✓ **Approved acc. to R129**



## **Trust in experience**

**20 years of experience** developing and manufacturing the best products.

- ✓ The only rear-facing rehab car seat with **full European approval**



## **Made in Germany**

All our products are developed and **manufactured in Germany.**

- ✓ **Certified according to ISO 9001**





## The big rear-facing car seat

The HERO-XL impresses with its ability to transport larger children in a rear-facing mode. This is usually the safest and most optimal position for children with limited trunk and head control.


The adapted padding system ensures also perfect support and comfort.

1	Total height	65 - 82 cm
2	Total width	43 cm
3	Total depth	55 cm

## As individual as your child

HERNIK offers various accessories to customise the HERO-XL even better to the needs of your child. For example, the car seat can be fitted with an **abduction block (4)**, a **seat depth adapter (5)** and a **security tray (6)**.

Children with weak trunk control usually have a need for a firm, stable support. The aluminium **trunk supports (7)** make this possible. They are completely flexible in width and height to the child's body. Unavoidable forces during acceleration, braking and cornering are absorbed so that your child feels safe and secure.



*We are very proud of the many positive customer ratings. On our website you can find out, why Hernik regularly achieves top ratings on comparison portals.*

**Markus Jacksch,**  
Managing Director



## Installation using a belt system

The HERO-XL does not require an ISOFIX connection and is secured in the vehicle using belts, which enables rearward travel up to 125 cm and 36 kg.

It can therefore be installed in almost any vehicle. The HERO-XL is the only rear-facing rehab car seat to date with full European approval in accordance with R129 (approval number: E4-08025600 S).



## Clear safety

### 5-point harness including padding ①

The integrated 5-point harness secures your child and ensures maximum safety in the event of an accident. It is very easy to operate using a quick-release fastener.

### Dynamic headrest technology ②

The height of the dynamic headrest can be infinitely adjusted with a quick-release fastener and thus individually adapted to the size of your child.

### Side impact protection ③

The side impact protection increases safety in the event of a side impact and is fitted on the side facing the vehicle door, leaving more space in the interior.

We are certain that we can only develop the best seats together with our customers. That is why we are very happy to receive your suggestions for improvement and will incorporate them into our development.

Please send us a message:  
[office@hernik.de](mailto:office@hernik.de)





## Further products from HERNIK




### Child car seats for every size

HERNIK accompanies your child from infancy to adulthood (up to 75 kg). For more complex requirements and individual requests we also offer customised solutions.

### Head protection helmets and seat shells

With our wide range of head protection helmets we offer every user the necessary protection against injuries. Our CAPOVARIO in particular is a real highlight that makes everyday life much easier thanks to its special design.

We also supply you with individual seat shells, such as the ANATOMIC UNICO and ANATOMIC FLEX aluminium shells.



# Hernik

**Discover all our products  
on [hernik.de](https://hernik.de)**

**HERNIK GmbH**

Gewerbering 8,  
83543 Rott am Inn

+49 (08039) 9093930  
[office@hernik.de](mailto:office@hernik.de)  
[www.hernik.de](https://www.hernik.de)

