


# IPA1 MAX

**Our Innovative Solution**

## Hernik



 ~ 4 - 12 years

100 - 150 cm



### **Comfortable & secure**

We rely on the **latest technology** to make our products safe.

✓ **Fast switch** between active and passive **seating position**



### **Trust in experience**

**20 years of experience** developing and manufacturing the best products.

✓ Rehab car seat with **full European approval**



### **Made in Germany**

All our products are developed and **manufactured in Germany**.

✓ **Certified according to ISO 9001**

## Adapts to the needs of your child!

The new IPAI-MAX is a customized car seat from Britax-Römer that can be tailored to meet the specific needs of your child and adapted accordingly.

Given that transferring often poses a significant challenge for parents, we recommend the entirely newly designed **rotating mechanism (1)**, which facilitates a 90° rotation. This innovation greatly simplifies the process of positioning the child in the car seat and securing the harness.

The torso guidance has also been taken into account: Aluminum **thorax pads (2)** are utilized, allowing for tool-free adjustments in both width and height—providing additional space in winter and a closer fit to the body in summer. An integrated **5-point harness (3)** is, of course, included as well.



Certainly, other prevalent accessories, including an **abduction block (4)**, a **foot rest (5)**, a **seat depth extension (6)**, and a **support and security tray (7)**, are also accessible.

The new IPAI-MAX has successfully undergone a comprehensive crash test in accordance with the R129 standard and possesses European approval with an E-number.

*We are very proud of the many positive customer ratings. On our website you can find out, why Hernik regularly achieves top ratings on comparison portals.*

**Markus Jacksch,**  
Managing Director

## Easily adjustable seating position

A completely redesigned ISOFIX **quick adjustment (1)** enables you to transition between an active and a passive seating position in mere seconds, as the **back angle (2)** is also infinitely adjustable. This can be accomplished swiftly and effortlessly without the need for any tools!

**This occurrence is unprecedented in the rehabilitation sector concerning car seats!**



Additionally, we are introducing a novel mechanism for **seat angle adjustment (3)** for the first time. Many children utilize a seat that is inclined backward in their wheelchair or buggy.

They can now adopt this familiar position in their car seat as well. With a single motion, the seat angle can be altered from flat to exceptionally steep—again, without the need for tools.

We are certain that we can only develop the best seats together with our customers. That is why we are very happy to receive your suggestions for improvement and will incorporate them into our development.

Please send us a message:  
[office@hernik.de](mailto:office@hernik.de)



## Further products from HERNIK




### Child car seats for every size

HERNIK accompanies your child from infancy to adulthood (up to 75 kg). For more complex requirements and individual requests we also offer customised solutions.

### Head protection helmets and seat shells

With our wide range of head protection helmets we offer every user the necessary protection against injuries. Our CAPOVARIO in particular is a real highlight that makes everyday life much easier thanks to its special design.

We also supply you with individual seat shells, such as the ANATOMIC UNICO and ANATOMIC FLEX aluminium shells.



# Hernik

**Discover all our products  
on [hernik.de](http://hernik.de)**

**HERNIK GmbH**

Gewerbering 8,  
83543 Rott am Inn

+49 (08039) 9093930  
[office@hernik.de](mailto:office@hernik.de)  
[www.hernik.de](http://www.hernik.de)

