

IPAI **MINI**

Our versatile companion

Hernik



~3 - 9 years

100 - 125 cm



Comfortable & secure

We rely on the **latest technology** to make our products safe.

✓ **Fast switch** between active and passive **seating position**



Trust in experience

20 years of experience developing and manufacturing the best products.

✓ Rehab car seat with **full European approval**



Made in Germany

All our products are developed and **manufactured in Germany**.


✓ **Certified according to ISO 9001**

Adapts to the needs of your child!

The new IPAI-MINI is a customized car seat from Britax-Römer that can be tailored to meet the specific needs of young children and adapted accordingly.

Given that transferring often poses a specific challenge for parents, we recommend the entirely newly designed **rotating mechanism (1)**, which facilitates a 90° rotation. This significantly simplifies the process of placing the child in the car seat and securing the harness.

The torso guidance has also been taken into account: **Aluminum thorax pads (2)** are utilized, which can be adjusted in width and height without the need for tools—providing more space in winter and a closer fit to the body in summer! An integrated **5-point harness (3)** is, of course, included.



Certainly, other prevalent accessories, including an **abduction block (4)**, a **foot rest (5)**, a **seat depth extension (6)**, a **lateral support (7)**, and a support and **safety cushion (8)**, are also accessible.

The new IPAI-MINI has successfully undergone a comprehensive crash test in accordance with the R129 standard and possesses European approval with an E-number.



We are very proud of the many positive customer ratings. On our website you can find out, why Hernik regularly achieves top ratings on comparison portals.

Markus Jacksch,
Managing Director

Easily adjustable seating position

A completely redesigned ISOFIX **quick adjustment (1)** enables you to transition between an active and a passive seating position in mere seconds, as the **back angle (2)** is also infinitely adjustable. This can be accomplished swiftly and effortlessly without the need for any tools!

This occurrence is unprecedented in the rehabilitation sector concerning car seats!



Additionally, we are introducing a novel mechanism for **seat angle adjustment (3)** for the first time. Many children utilize a seat that is reclined in their wheelchair or buggy. They can now adopt this familiar position in their car seat as well. With a single motion, the seat angle can be adjusted from flat to significantly steep—again, without the need for tools.

We are certain that we can only develop the best seats together with our customers. That is why we are very happy to receive your suggestions for improvement and will incorporate them into our development.

Please send us a message:
office@hernik.de



Further products from HERNIK



Child car seats for every size

HERNIK accompanies your child from infancy to adulthood (up to 75 kg). For more complex requirements and individual requests we also offer customised solutions.

Head protection helmets and seat shells

With our wide range of head protection helmets we offer every user the necessary protection against injuries. Our CAPOVARIO in particular is a real highlight that makes everyday life much easier thanks to its special design.

We also supply you with individual seat shells, such as the ANATOMIC UNICO and ANATOMIC FLEX aluminium shells.



Hernik

**Discover all our products
on hernik.de**

HERNIK GmbH

Gewerbering 8,
83543 Rott am Inn

+49 (08039) 9093930
office@hernik.de
www.hernik.de

