

Am Eckfeld II, Gewerbering 8
83543 Rott am Inn
Germany
Tel.: +49 (0)8035-9580448
Fax: +49 (0)8035-9580446
info@hernik.de • www.hernik.de



invoice / delivery address

Kom.:

☐ **Order**

☐ **Offer**

Rehab Car Seat

KIDSFLEX

R129

NEW!

Edition

Exportprices

Release: 2025

BASIC

	size 1	size 2	size 2XL
total height	74 cm	80 cm	90 cm
crown height	67 cm	73 cm	82 cm
back height	37 - 51 cm	41 - 57 cm	47 - 65 cm
seat width	31 cm	35 cm	40 cm
seat depth (adjustable)	28 - 33 cm	34 - 39 cm	40 - 46 cm
lateral support height (small)	26 - 36 cm	26 - 36 cm	34 - 44 cm
lateral support height (large)	30 - 40 cm	30 - 40 cm	38 - 48 cm
distance of lateral supports	16 - 28 cm	19 - 32 cm	24 - 37 cm
weight	8,2 kg	8,9 kg	12,3 kg
capacity	50 kg	50 kg	100 kg

common information (standard):

NEW!

- seat depth adjustable
- headrest included (longer backrest)
- upholstery & cover
- backrest angle adjustable
- 5-point-harness with upholstery
- frame to fix a footrest, clamp bow or ISOFIX
- Crash tested according to **R129/03**
- **ISOFIX adapter** included
- **headrest KICK-BACK** included
- **trunk supports** small, fix, included
- **shoulder support** included



- ☐ size 1
black
(300000-E3)
- ☐ size 2
black
(310000-E3)
- ☐ size 2XL
black
(310000-E3-XL)

belt guide*:

- ☐ LEFT
- ☐ RIGHT

*please choose ONE of them

- ☐ wider hip
cushions (SW -4 cm)
(300850-E)

fittings for the different needs

head support „Kick-Back“
with aluminum frame
(free of charge)



☐ head fixation
(300630-E)



leg support

☐ left
(300350-E)

☐ right
(300355-E)



abduction block, removable

☐ size 1
(300400-E)

☐ size 2
(310400-E)

☐ size 2XL
(310410-E)




foot rest, height and angle adjustable

☐ size 1
(300450-E)

☐ size 2, 2XL
(310450-E)


☐ swing up
(300452-E)



foot roll, height and angle adjustable

☐ size 1
(300455-E)

☐ size 2, 2XL
(310455-E)



trunk supports, small
(height: 11 cm, depth: 14 cm)

☐ fix
(free of charge)

☐ swing away
(300505-E3)



trunk supports, large
(height: 21 cm, depth: 16 cm)

☐ fix
(300510-E3)

☐ swing away
(300515-E3)



shoulder supports
(free of charge)



 fittings frame

swivel base

☐ size 1 and 2
(300200-E)

☐ RIGHT

☐ LEFT



fastener

(independent of the direction of rotation)

☐ inside
(standard)


☐ outside
(300220-E)

angular adjustment (0° - 10°)

☐ size 1
(300100-E)

☐ size 2
(310100-E)

☐ size 2XL
(310100-E-XL)



extra angular adjustment (10° - 20°)

☐ size 1
(300105-E)

☐ size 2
(310105-E)

☐ size 2XL
(310105-E-XL)



☐ Sandwich-Base (mounted under the base board) (300790-E)

Useful in the following situations:

- with very sporty shaped car seats, in order to compensate for the recess there. This results in a more stable position and more freedom of movement when using the swivel adapter.
- with car Isofix brackets that are very far inside to allow nevertheless contact with the car's backrest and the KIDSFLEX (the ISOFIX adapter sits in the sandwich base)
- approx. 1.5 kg additional weight


fittings for car installation

Fixation (one option included in delivery)

ISOFIX-Adapter

 size 1
(300575-E3)


 size 2 & 2XL
(310575-E3)



clamp bow


 size 1
(300550-E3)

 size 2 & 2XL
(310550-E3)



outrigger*

☐ (310453-E)



*300450-E
or 310450-E
necessary!

additional belt guide
(on the other side)

☐ (300075-E)



fittings security

security tray

 size 1
(300600-E)

 size 2
(310600-E)

 size 2XL
(310600-E)



security pad for the buckle
with snap fastener

 size 1, 2, 2XL
(300650-E)

The buckle gets invisible for the child thanks to an additional pocket!



security pad for the buckle
with magnetic lock

 size 1, 2, 2XL
(300655-E)



protection mat
(to protect the original seats of the car)

☐ (300780-E)



individual colors for the cushions
(300800-E)

[one X means single-colored, two X means bicolored]

	inside	outside
☐ red	☐	☐
☐ blue	☐	☐
☐ black	☐	☐
☐ grey	☐	☐

belt to connect the chest straps
(to improve the position of the belt in front of the breast)

☐ (300640-E)



notes