

Am Eckfeld II, Gewerbering 8
83543 Rott am Inn
Tel.: +49 (0)8039-9093930
Fax: +49 (0)8039-90939319
info@hernik.de • www.hernik.de

invoice / delivery address

Kom.:

☐ **Order**
☐ **Offer**

Rehab Car Seat **KIDSFLEX^{PRO}**

Exportprices

Release: 2026

BASIC

	size 1	size 2
crown height	56 - 71 cm	58 - 73 cm
back height	33 - 48 cm	33 - 51 cm
seath width	28 cm	33 cm
seat depth (adjustable)	31 - 37 cm	36 - 42 cm
height trunk supports, small	19 - 38 cm	24 - 43 cm
height trunk supports, large	23 - 42 cm	28 - 47 cm
distance trunk supports	14 - 25 cm	17 - 28 cm
body height	100 - 125 cm (max. 50 kg)	100 - 150 cm (max. 50 kg)
total height	71 - 86 cm	72 - 87 cm



common information (standard):

- crash tested according to ECE-R129/04
- 5-point-harness with upholstery, height adjustable
- head support „standard“ inclusive
- base board with ISOFIX adapter inclusive
- stepless back angle adjustment
- adjustable seat depth
- shoulder supports inclusive
- trunk supports, small inclusive
- adjustable back plate


☐ size 1
grey-black
(360000-E1)

☐ size 1
red-black
(360000-E2)

☐ size 2
grey-black
(370000-E1)

☐ size 2
red-black
(370000-E2)

belt guide (please ALWAYS choose one of them)

☐ LEFT
(300000-L)

☐ RIGHT
(300000-R)

☐ seat cushion
standard

(350355-E)

☐ seat cushion PRO (ergonomically shaped
seat cushion including abduction pads and
width-adjustable leg supports)

 size 1 (360350-E)
size 2 (370350-E)

fittings for different needs

head support, standard
(free of charge)



☐ head support PRO
(350025-E)



☐ head fixation
(350630-E)



☐ wider hip cushions
(reduces the seat width by 4 cm)

(350850-E)

abduction block, removable

☐ size 1
(360410-E)

☐ size 2
(370410-E)



foot rest, height and angle adjustable


☐ size 1 and 2
(350450-E)

☐ swing up
(350452-E)



foot roll, height and angle adjustable


☐ size 1 and 2
(350455-E)



trunk supports, small
(height: 11 cm, depth: 15 cm)

☒ fix
(350500-E)

☐ swing away
(350505-E)



trunk supports, large
(height: 21 cm, depth: 16 cm)

☐ fix
(350510-E)

☐ swing away
(350515-E)



shoulder supports
(free of charge)


 fittings frame

swivel base

☐ (350200-E)

☐ RIGHT

☐ LEFT



☐ angular adjustment (5° - 20°)

size 1 (360100-E)

size 2 (370100-E)



☐ conversion set

☐ from ISOFIX adapter
to clamp bow
(350020-E)

☐ from clamp bow
to ISOFIX adapter
(350030-E)

fittings for car installation

Fixation (one option included in delivery)

ISOFIX adapter

☐ size 1
(360575-E)


☐ size 2
(370575-E)


clamp bow

☐ size 1
(360550-E)

☐ size 2
(370550-E)


outrigger*

☐ (350453-E)

*350450-E
or
350450-E
necessary!

additional belt guide
(on the other side)

☐ (350075-E)

 fittings security

security tray
(in height and depth adjustable)

☐ size 1
(360600-E)

☐ size 2
(370600-E)

security pad for the buckle
Cover with snap fastener

☐ (350650-E)

security pad for the buckle
with HERNIK Security Key*

☐ (350655-E)

NEW!

* incl. belt cutter

protection mat
(to protect the original seats of the car)

☐ (350780-E)

individual colors for the cushions
(300800-E)

(one X means single-colored, two X means bicolored)

	inside	outside
<input type="checkbox"/> red	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> black	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> grey	<input type="checkbox"/>	<input type="checkbox"/>

belt to connect the chest straps

☐ (350640-E)

 notes
